LSSTH BHAVANAGAR +91-278-2423250

EQAS CHECKLIST

UNIQUE ID LSSTH/ BIOCHEM/EQASCHECKLIST PAGE NO 1 of 3

Danard	
Department: Biochemistory	EQAS Program: RICAS
10.	Cycle: 18
17201100.	Sample: 10
Lab Result: 89.20	Instrument: AFT 300
Z-score:	Peer Mean:
Comparator Mean: 94.86	Comparator (n):

Sample Receipt:	Who Received Kit:	Date Received:	Yes/No/NA	
Sample Receipt:	Pallar	25/03/2021	e Romania di	
Was the kit received in	good condition?	The state of the s	428	
Was the Kill received in	good condition:		Yes	
Was the correct program/cycle received? Was the kit stored at the proper temperature following receipt?			yes.	
Sample Preparation:	Who Prepared Sample:	Date: 22/10/2021		
	itution instruction followed (h		res	
temperature)?	A (or calibrated) pipette used	for reconstitution?	ves	
was a volumetric Class	zed water used for reconstitution	on?	yes yes yes	
Was the sample mixed	according to the package inser	rt prior to testing?	yes	
Sample Processing:	Who Ran Affected Test:	Date: 22/10/2021		
Was the correct sample number tested?			Yes	
- Was the sample at roo	om temperature?		Yes	
- Was the person running the test current in their training			No	
 Was the affected test 	run within the stability claim	listed in the package	Yes	
insert? Was incubation time for analysis is followed correctly?			yes	
Reporting Results:	Who Reported?	Date: 25/10/2021		
- Was the test configuration correct (instrument, method and reagent)?			Yes	
- Have the results been reported correctly (match instrument print out)?		74		
- Was the correct unit reported?		Yes		
- Was the decimal symbol placed correctly when reported?			1es	

Name of Laboratory : Labo	ratory Services Sir T. Hospital		
Document No.5	Document Name : SOP C	F IQC and EQAS	
Issue No. :01	Issue Date: 26/10/2015	Amend No.:	Amend Date :
Prepared by: Section In charge		Approved b HOD,Bioch	4

LSSTH BHAVANAGAR

+91-278-2423250

EQAS CHECKLIST

UNIQUE ID LSSTH/ BIOCHEM/EQASCHECKLIST PAGE NO 2 of 3

- Was the reported resu	It within the instrument's linea	r range?	Yes
- Was the calculation of	f the reported result done corre	ctly?	
Internal QC:	Who Ran IQC (on EQAS testing):	Date: 07	. 185
Was the correct reconst	itution instruction followed (ha		
temperature)?			Yes
Was IQC material store	d at the proper temperature fol	lowing receipt?	yes
Was IQC within an acc	eptable range on the day that th	ne EQAS sample was run?	YRS
Where there any shifts run?	or trends in IQC just before/aft	er the EQAS sample was	No
Calibration:	Who Ran Calibration	Date: 92/10/2021	
Was the last calibration	acceptable?	*	yes
Was the last calibration	within the manufacturer's reco	ommended dating?	Yes
Reagent:	Reagent Lot: DE 2201001	Lot Expiration 28/06/2022	
-Was the test reagent s			Yes
- Was the test reagent properly prepared?			Yes
- Was the test reagent within manufacturer's dating?			yes
Instrument:	Who Performed Maintenance: Leding	Date: 22/10/202)	
Was daily maintenance	e performed on the day that the		Yes
	ning maintenance current on tr		No
Was the instrument op	erating correctly on the day the	sample was tested?	y.e3
	ent acceptable for the instrumer		yes
Was instrument calibra	ated or not?		
Sample Retest:	Who Retested Sample:	Date: 7/11/2021	
Was the EQAS sample	retested following receipt of E	EQAS sample report?	Yes
If yes, was the result within acceptable limits for the EQAS sample?			Yes C92-26
EQAS Evaluation:			
Inappropriate peer gro	up or comparator	8 8 2	
Inappropriate evaluation method or instrument.	on criteria, e.g. narrow limits d	ue to use of a precise	
Name of Laboratory · Labo	ratory Services Sir T. Hospital		
Document No.5	Document Name : SOP OF IQ	OC and EOAS	and the state of the state of
Issue No. :01			
Prepared by: Section In charge		Approved by: HOD,Biochemistry	nd Date :

LSSTH DIV		TO LCCTU/
	EQAS CHECKLIST	UNIQUE ID LSSTH/ BIOCHEM/EQASCHECKLIST
+91-278-2423250		PAGE NO 3 of 3

- -> After receiving the EOAS result, we had repeated (Refested) the ECLAS Semple for serum chloside.
 - > We found that probable course of outlier may be due to Remdom orran.

Corrective Actions:

-> Pe tested the ECRAS Sumple and Fresult found Sectisfactory.

Preventive Actions:

We will check for any Irundom error before Trynning the ECLAS and other samples.

Reviewed By:

Japan Dr. Supra Partel

Name of Laboratory . Dabor	atory Services Sir T. Hospital			St. A. Carlotte and Carlotte	A Part of the second
Document No.5	Document Name: SOP OF IQC and EQAS				
Issue No.:01	Issue Date: 26/10/2015	Amend No.:	_	Amend Date:	
Prepared by: Section In charge	2		proved by: OD,Biochemistr	ry 🚮	

